




Supplement of

A novel inlet system for online chemical analysis of semi-volatile submicron particulate matter

P. Eichler et al.

Correspondence to: A. Wisthaler (armin.wisthaler@uibk.ac.at)

SUPPLEMENTARY INFORMATION:



Valve settings:

	V1	V2	V3	V4	V5
particle phase	open	closed	open	closed	closed
particle phase (zero)	open	closed	closed	open	open
gas phase	closed	open	closed	open	open
gas phase (zero)	gas phase inlet is overflowed with catalytically cleaned sampling air				