



*Supplement of*

## **Spectral Aerosol Extinction (SpEx): a new instrument for in situ ambient aerosol extinction measurements across the UV/visible wavelength range**

C. E. Jordan et al.

*Correspondence to:* C. E. Jordan ([carolyn.jordan@nianet.org](mailto:carolyn.jordan@nianet.org))

The copyright of individual parts of the supplement might differ from the CC-BY 3.0 licence.

## **SUPPLEMENTARY TABLES**

**Table S1.** Minimum standard deviations for select wavelengths from Allan Variance test. Test conducted with filtered air at a flow rate of 90 lpm using 50 ms integration time per spectrum.

<i>Wavelength (nm)</i>	<i># of scans for minimum st. dev.</i>	<i>Time (s) to acquire # of scans</i>	<i>Minimum st. dev. (Mm<sup>-1</sup>)</i>
300	1379	69	5.7
400	1076	54	3.5
500	946	47	4.2
600	946	47	5.4
700	946	47	3.9

**Table S2.** Summary of aerosol tests: materials used and generated size distributions. A total of 22 tests were performed generating 87 total spectra.

<i>Material</i>	<i>Size Distributions</i>			
	<i>200 nm</i>	<i>600 nm</i>	<i>900 nm</i>	
<i>Scatterers</i>				
Polystyrene Latex Spheres	2	1	1	
Ammonium Sulfate	1	1		
<i>Dust</i>				
Montmorillonite	1	2		
Blue Ridge Violet	1			
Italian Yellow Earth	1			
Luberon Natural Red	1	1		
Goethite	1			
<i>Smoke &amp; Soot Analogs</i>				
Fullerene Soot				
Aquadag				
Smoke		1	1	
Low O:C Smoke	1			
High O:C Smoke	1			

## **SUPPLEMENTARY FIGURE CAPTIONS**

**Figure S1.** Intensity spectrum of the xenon lamp over the full wavelength range of the spectrometer for a filtered air sample. The measurement range is highlighted by the red box.

**Figure S2.** Sensitivity of MiePlot  $C_{ext}$  calculations to particle diameter for monodisperse PSL spheres (top panel) and to lognormally distributed PSL spheres with standard deviations of 1% and 2% (bottom panel).




Figure S1.





Figure S2.