

1 **Development of an SI-traceable transmission curve reference material** 2 **to improve quantitation and comparability of PTR-MS measurements**

3 David R. Worton¹, Sergi Moreno^{1,#}, Kieran O'Daly¹, Rupert Holzinger²

4 ¹ National Physical Laboratory, Hampton Road, Teddington, TW11 0LW, United Kingdom

5 ² Institute for Marine and Atmospheric Research, IMAU, Utrecht University, Utrecht, the Netherlands

6 *Correspondence to:* David R. Worton (dave.worton@npl.co.uk)

7 [#] Now at World Meteorological Organization (WMO), 7bis Avenue de la Paix, C.P. 2300

8 CH-1211, Geneva 2, Switzerland

9

10

11 **Abstract.** Since its inception more than two decades ago proton-transfer-reaction mass-spectrometry (PTR-MS) has
12 established itself as a powerful technique for the measurements of a wide range of volatile organic compounds (VOCs) with
13 high time resolution and low detection limits and without the need for any sample pre-treatment. As this technology has
14 matured and its application become more widespread there is a growing need for accurate and traceable calibration to ensure
15 measurement comparability. As a result of the large number of VOCs detectable with PTR-MS it is impractical to have a
16 calibration standard or standards that cover all observable compounds. However, recent work has demonstrated that
17 quantitative measurements of uncalibrated compounds are possible provided that the transmission curve is accurately
18 constrained. To enable this, a novel traceable multi-component gas reference material containing 20 compounds spanning a
19 mass range of 32 to 671 has been developed. The development and compositional evolution of this reference material is
20 described along with an evaluation of its accuracy and stability. This work demonstrates that for the majority of components
21 the accuracy is < 5 % (most < 3 %; < 10 % for hexamethylcyclotrisiloxane (D3-siloxane) and 1,2,4-trichlorobenzene (1,2,4-
22 TCB)) with stabilities of > 2 years (> 1 year for acetonitrile, methanol and perfluorotributylamine (PFTBA)).

23

24 **1 Introduction**

25 Proton-transfer-reaction mass-spectrometry (PTR-MS) is a technique that allows simultaneous measurements of multiple
26 volatile organic compounds (VOCs) in real-time (≤ 1 s) with low detection limits (pmol mol^{-1}) and without any sample pre-
27 treatment (Lindinger et al., 1998; Hansel et al., 1999). Most VOCs, with the exception of alkanes with less than five carbon
28 atoms, have proton affinities larger than water ($691 \pm 3 \text{ kJ mol}^{-1}$) and are therefore detectable with PTR-MS. Also, because
29 most VOCs have proton affinities below 900 kJ mol^{-1} there is minimal excess energy following proton transfer resulting in
30 minimal fragmentation. For these reasons it is a very convenient measurement technique for a wide range of applications. Over
31 the last two decades PTR-MS has become an important and widely applied tool for VOC measurements that has resulted in
32 major advances in the field of atmospheric sciences (De Gouw and Warneke, 2007; Park et al., 2013; Yuan et al., 2017). It has
33 also been applied in the medical sector for the detection of VOCs to diagnose diseases or disease states (Beauchamp et al.,
34 2013) and in the food and beverage industry for characterising flavour and odour (Biasioli et al., 2011).

35

36 Multiple manufacturers now produce and commercially sell PTR-MS instruments globally that differ in the production and
37 detection of ions including different types of mass spectrometer. Therefore, as its application becomes more widespread, and
38 more datasets are generated there is an increasing need for accurate calibration and measurement comparability. Additionally,
39 as part of the European funded Aerosols, Clouds and Trace Gases Infrastructure (ACTRIS) project (<https://www.actris.eu/>)
40 there is an interest to establish PTR-MS as a technique for long-term monitoring of VOCs, which emphasises the need for a
41 robust metrological infrastructure to control and assure the quality of data produced by monitoring stations performing these
42 measurements. However, the lack of traceable reference materials to calibrate PTR-MS instruments presents challenges in the
43 pursuit of obtaining comparable results and is an obstacle to long-term studies. Primary reference materials (PRMs) prepared
44 by gravimetry in high pressure cylinders by national metrology institutes (NMIs) underpin the accuracy (trueness) and
45 comparability of measurement data through traceability to the international system of units (SI). Traceability has been
46 demonstrated as a critical component for chemical measurements that ensures the comparability, stability and coherence in
47 measurements providing confidence in measurement results (Brown and Milton, 2007). PRMs produced by NMIs represent
48 the highest point in the traceability chain and the accuracy and international comparability is ensured through key comparisons
49 organised within the Consultative Committee on Amount of Substance Gas Analysis Working Group (CCQM-GAWG) and in
50 regional comparisons organised within the Regional Metrology Organisations (RMOs), e.g., EURAMET (Europe).

51

52 As a result of the numerous VOCs detectable with PTR-MS it is impractical to have a calibration standard or standards that
53 cover all observable compounds. However, since the conception of PTR-MS, there has been awareness for the potential of this
54 technique to provide quantitative measurements for compounds without the need for specific calibration materials (Hansel et
55 al., 1999). The basis for this is that the amount fraction of compound R ($[R]$) can be determined from (Taipale et al., 2008):

56

$$[R] = \frac{1}{k\Delta t} \times \frac{I(RH^+)}{T(RH^+)} \times \left(\frac{I(H_3O^+)}{T(H_3O^+)} \right)^{-1} \quad \text{Eq. 1}$$

58

59 Where k is the proton transfer reaction rate coefficient, Δt is the reaction time, $I(RH^+)$ and $I(H_3O^+)$ are the observed ion count
60 rates for the protonated ion of compound R (RH^+) and the hydronium ion (H_3O^+), respectively. $T(RH^+)$ and $T(H_3O^+)$ are the
61 transmission efficiencies for RH^+ and H_3O^+ ions, respectively. The transmission coefficients are predominantly mass
62 dependent, but they can also vary in time (De Gouw et al., 2003; Ammann et al., 2004; Steinbacher et al., 2004). Proton transfer
63 reaction rate coefficients can be measured and/or predicted using quantum methods (Zhao and Zhang, 2004). If specific rate
64 coefficients are agreed within the community for specific compounds and are widely used this would negate the role of different
65 rate constants on measurement comparability (Table S1, Supporting Information). The reaction time and observed ion count
66 rates are all measured parameters leaving just the transmission coefficients as variables required for quantitative measurements
67 without specific calibrations. Cappellin et al. (2012) demonstrated the quantitative properties of one type of PTR-MS
68 instrument by assuming a theoretical transmission based on the duty cycle of the time-of-flight mass analyser. However, for
69 newer generation instruments that employ advanced ion optics to improve sensitivity, it is necessary to determine the mass-
70 dependent transmission experimentally as the transmission of the system diverges from theory at low masses. Deviations can
71 also occur at high masses due to poor tuning and/or ageing of the ion detection system (Müller et al., 2014).

72

73 There are several highly cited publications that explore best practices in PTR-MS measurements (e.g., Blake et al., 2009; De
74 Gouw and Warneke, 2007; Yuan et al., 2017), including methods to calibrate and retrieve the mass dependent transmission
75 (Taipale et al., 2008). However, many of these methods are slow and labour intensive and as a result calibrations and
76 transmission curve retrievals are not performed frequently enough. This has limited the application of PTR-MS to mostly short
77 campaign-scale intensive deployments and only a few groups have utilised PTR-MS for long-term studies (Holzinger et al.,
78 2006; Taipale et al., 2008). However, recent work by Holzinger et al., (2019) has demonstrated: (i) a new method to retrieve
79 the mass-dependent transmission from fast calibrations that should enable more frequent calibrations and (ii) the validity of a
80 simple reaction kinetics approach to quantify measurements of uncalibrated compounds from different PTR-MS instruments
81 with an accuracy of $\leq 30\%$ provided the transmission curve is accurately constrained. A prototype PRM, 0917a reported in
82 this work, was initially developed by the National Physical Laboratory (NPL) and employed for the PTR-MS intercomparison
83 campaign at the CESAR observatory in the central Netherlands (Holzinger et al., 2019). Following this comparison exercise
84 improvements to the composition were needed to include additional compounds in the mass-to-charge (m/Q) 150 – 400 Th
85 range to provide a more robust retrieval of the mass-dependent transmission.

86

87 In this paper, the development and compositional evolution of PRMs and certified reference materials (CRMs) specific for
88 constraining the PTR-MS transmission curve are described, including an evaluation of the accuracy through comparisons

89 validating the gravimetric preparation of various different PRMs of similar composition and an assessment of their long-term
90 stabilities. For details on how to use the RMs to constrain the PTR-MS transmission curve the reader is directed to H

91 **2 Experimental methods**

92 **2.1 Gravimetric preparation of primary reference materials**

93 The PRMs were prepared at four distinct timepoints (September 2017, December-January 2018, August 2019 and August
94 2021) and the compositions evolved over this timeframe (Table 1) due to improvements in the preparation and validation
95 techniques (e.g., 1,2,4-trichlorobenzene; 1,2,4-TCB) or due to requests from the PTRMS community for inclusion of new
96 components (e.g., dimethyl sulfide; DMS). All the PRMs were prepared gravimetrically in accordance with ISO 6142-1:2015
97 (ISO, 2015) from pure components. All pure components were purity analysed in accordance with ISO 19229 (ISO, 2019).
98 Table S2 (Supporting Information) provides the sources and purities for each component and shows that all chemicals with
99 the exception of perfluorotributylamine (PFTBA) were $\geq 98\%$ pure. Table S2 (Supporting Information) also shows the boiling
100 points and vapour pressures for all compounds. All components were liquids at room temperature and pressure, with the
101 exception of propane (gas) and hexamethylcyclotrisiloxane (D3-siloxane; solid). As a solid under room temperature and
102 pressure conditions, D3-siloxane needed to be dissolved into a solvent to enable its addition to the cylinder. Further details are
103 given in the Supporting Information (Supplementary text: Preparation of D3-siloxane reference materials).

104

105 All PRMs were prepared in 10 L aluminium cylinders (Luxfer) with a proprietary passivation treatment (Spectraseal™, BOC)
106 and BS341 no. 15 outlet diaphragm valves (Ceodeux). Cylinders were evacuated using an oil free pump (Scrollvac SC15D,
107 Leybold Vacuum) and turbo molecular pump with magnetic bearing (Turbovac 340M, Leybold Vacuum) to a pressure of < 3
108 $\times 10^{-7}$ mbar. Individual compounds were added to the evacuated cylinder via a transfer vessel (capped $1/8$ " diameter tube, with
109 a nominal volume of 1 mL, Swagelok, electro-polished stainless steel). The transfer vessel was weighed against a tare vessel
110 matched for size and shape before and after each addition into the evacuated cylinder (Mettler-Toledo XP2004S). The ultra-
111 high purity nitrogen balance gas (BIP⁺, Air Products) was added via direct addition to the cylinder, through purged $1/16$ " tubing
112 (Swagelok, electro-polished stainless steel). For the vast majority of compounds, they were initially produced as binary parent
113 mixtures at amount fractions $> 10 \mu\text{mol mol}^{-1}$ (typically at nominally $50 \mu\text{mol mol}^{-1}$) though some were produced as ternary
114 or quaternary mixtures containing two or three compounds together in the same parent mixture. A full breakdown of the 50
115 parent mixtures used to prepare the six PRMs developed in this work are shown in Table 2. Aliquots of each of these parent
116 mixtures were added by direct addition to an evacuated cylinder to produce a final mixture containing all 16 to 20 VOCs at
117 nominal amount fractions of $1 \mu\text{mol mol}^{-1}$ in a balance of nitrogen.

Table 1. Overview of composition (name, formula, CAS#), the protonated monoisotopic molecular ion and any major fragment ions, formed following protonation in H₃O⁺ mode) for 20 compounds included in the PRM and CRM reference materials and the key comparisons through which the traceability to the international community is derived.

Compound ^a	Formula	CAS #	m/Q [Th]		Traceability (Reference, if applicable)
			protonated	fragments	
methanol	CH ₃ OH	67-56-1	33.033	-	EURAMET-1305 OVOCs (Worton et al., 2022), CCQM-K174 ^b
acetonitrile	CH ₃ CN	75-05-8	42.034	-	CCQM-K174 ^b
acetaldehyde	CH ₃ CHO	75-07-0	45.033	-	CCQM-K174 ^b
propane ^c	C ₃ H ₈	74-98-6	not detected	-	CCQM-K111 (Veen et al., 2017)
ethanol	C ₂ H ₅ OH	64-17-5	47.049	-	CCQM-K93 (Brown et al., 2013), EURAMET-1305 OVOCs (Worton et al., 2022), CCQM-K174 ^b
acetone	CH ₃ OCH ₃	67-64-1	59.049	-	EURAMET-1305 OVOCs (Worton et al., 2022), CCQM-K174 ^b
DMS	C ₂ H ₆ S	75-18-3	63.026	-	CCQM-K94 (Lee et al., 2016), CCQM-K165 (Lee et al., 2022)
isoprene	C ₅ H ₈	78-79-5	69.070	41.039	EURAMET-886 (Grenfell et al., 2008; Grenfell et al., 2010)
MVK	C ₄ H ₆ O	78-94-4	71.049	-	CCQM-K174 ^b
MEK	C ₄ H ₈ O	78-93-3	73.065	-	CCQM-K174 ^b
benzene	C ₆ H ₆	71-43-2	79.054	-	CCQM-K10.2018 (Cecelski et al., 2022)
toluene	C ₇ H ₈	108-88-3	93.07	-	CCQM-K10.2018 (Cecelski et al., 2022)
m-xylene	C ₈ H ₁₀	108-38-3	107.086	-	CCQM-K10.2018 (Cecelski et al., 2022)
1,2,4-TMB	C ₉ H ₁₂	95-63-6	121.101	-	EURAMET-886 (Grenfell et al., 2008; Grenfell et al., 2010)
1,2,4-TFB	C ₈ H ₈ F ₃	367-23-7	133.026	-	-
3-carene	C ₁₀ H ₁₆	13466-78-9	137.132	81.070	CCQM-K121 (Liaskos et al., 2018)
1,2,4-TCB	C ₆ H ₅ Cl ₃	120-82-1	180.937	-	-
D3-siloxane ^d	C ₆ H ₁₈ O ₃ Si ₃	541-05-9	223.064	207.033, 225.033 ^e	EURAMET 1305 Siloxanes ^b
D4-siloxane ^e	C ₈ H ₂₄ O ₄ Si ₄	556-67-2	297.082	281.051, 283.030, 299.062 ^e	EURAMET 1305 Siloxanes ^b
D5-siloxane ^f	C ₁₀ H ₃₀ O ₅ Si ₅	541-02-6	371.101	355.070, 373.081 ^e	EURAMET 1305 Siloxanes ^b
PFTBA	C ₁₂ F ₂₇ N	311-89-7	not detected	651.961, 413.977	-

^aShort names shown here but preferred IUPAC names are: propan-2-one (acetone), (Methylsulfonyl)methane (dimethyl sulfide; DMS), 2-Methylbuta-1,3-diene (isoprene), but-3-en-2-one (methyl vinyl ketone; MVK), butan-2-one (methyl ethyl ketone; MEK), 1,2,4-trimethyl benzene (1,2,4-TMB), 1,2,4-trifluoro benzene (1,2,4-TFB), 3,7,7-trimethylbicyclo[4.1.0]hept-3-ene (3-carene), 1,2,4-trichlorobenzene (1,2,4-TCB), hexamethylcyclotrisiloxane (D3-siloxane), octamethylcyclotetrasiloxane (D4-siloxane), decamethylcyclopentasiloxane (D5-siloxane), and perfluorotributylamine (PFTBA). ^bComparison in progress at the time of publication. ^cNot detectable in PTRMS in H₃O⁺ mode but included as an internal standard. ^dFurther information on the mechanisms yielding product ions and fragments in Fig. S4 (Supporting Information). ^eMethyl/hydroxyl group switch.

Table 2. Composition, hierarchies and parent cylinder IDs (dates prepared) for all components for the six PRMs (0917a, 0917b, 1218, 0119, 0819 and 0821) prepared in this work. The PRMs are identified by the date and year of their preparation (MMYY). As the first two were produced at the same time the suffices a and b are added to distinguish them. The colour scheme in the table headers is matched to that used in the figures throughout the paper.

Compound	Cylinder ID (Date Prepared)					
	0917a (18/09/2017)	0917b (18/09/2017)	1218 (04/12/2018)	0119 (02/01/2019)	0819 (23/08/2019)	0821 (31/08/2021)
methanol		A463 (13/02/2015)			A410 (24/01/2013)	A602 (26/04/2018)
acetonitrile		A389R (11/02/2015)			A403 (11/10/2012)	A670R (01/04/2021)
acetaldehyde		A400R (02/02/2015)			A402 (11/10/2012)	2832 (01/04/2021)
propane		D91038IR (18/11/2014)			NG561 (16/10/2014)	-
ethanol		A463 (13/02/2015)			A410 (24/01/2013)	A602 (26/04/2018)
acetone		A463 (13/02/2015)			VOC6 (05/05/2009)	A602 (26/04/2018)
DMS	-		2106 (21/02/2017)		NG388 (13/09/2012)	3073 (16/11/2020)
isoprene		D292194R (13/01/2011)			VOC6 (05/05/2009)	D994138R2 (28/09/2020)
MVK		2064 (24/06/2016)			2088 (24/06/2016)	3070 (02/08/2021)
MEK		A389R (11/02/2015)			A403 (11/10/2012)	3070 (02/08/2021)
benzene		D91038IR (18/11/2014)			D842635R (13/10/2015)	D618317 (15/08/2018)
toluene						D600070 (19/03/2018)
m-xylene	D641688 (04/03/2010)		D618307 (19/03/2018)		D994171 (20/03/2013)	D618307 (19/03/2018)
1,2,4-TMB	D442684 (02/03/2017)			D711530 (26/11/2018)		D442684 (02/03/2017)
1,2,4-TFB		A569 (11/09/2017)			2810 (21/08/2019)	D723197R (14/06/2021)
3-carene	D090493 (18/11/2014)				D711532 (21/11/2018)	
1,2,4-TCB	-		D641970R (31/07/2018)		A568 (18/08/2017)	D723197R (14/06/2021)
D3-siloxane	-		2586 (07/11/2018)	2693 (07/11/2018)	2586 (07/11/2018)	3134 (22/06/2021)
D4-siloxane	A582 (15/08/2017)	A567 (18/08/2017)			A582 (15/08/2017)	A629R (22/06/2021)
D5-siloxane	A582 (15/08/2017)	A567 (18/08/2017)	A629 (26/11/2018)		A582 (15/08/2017)	A629R (22/06/2021)
PFTBA						D961497 (30/07/2021)

120 2.2 Analytical methods

121 To perform the analytical validation a method was developed on a gas chromatograph (GC; Agilent Technologies 7890)
122 instrument equipped with both a flame ionisation detector (FID) and electron ionisation (70 eV) mass spectrometer (Agilent
123 5975; GC-FID/MS system). Samples were introduced using a 6 port 2 position valve (VICI) and a fixed sampling loop (1 mL).
124 The column effluent was split to both detectors simultaneously by using a detector splitter plate (Agilent Technologies).
125 Separation was achieved for all components using a DB-624 capillary column (J&W; 75m × 0.53 mm, df = 3 µm). The carrier
126 gas was helium and the flow was held constant at 5 mL min⁻¹, with a temperature program starting at 30 °C held for 10 minutes,
127 ramped at 10 °C min⁻¹ to 120 °C and held for 15 minutes before a final ramp of 50 °C min⁻¹ was applied to a final temperature
128 of 200 °C, which was held for a further 10 minutes. The total run time was 46 minutes.


129

130 Low FID responses for methanol and acetaldehyde presented some analytical challenges because the observed peaks were too
131 small to achieve useable results due to poor reproducibility. As a result, another analytical method was developed on a second
132 GC-FID instrument without an MS (Scion 456; Cryo-GC-FID system) that had a pre-concentration trap (15 cm of 1/8" tubing;
133 1 mL volume) packed with glass beads and cooled with liquid nitrogen that enabled trapping of larger volume samples yielding
134 larger peaks and improved repeatability for all three compounds. The pre-concentration trap was held at -185 °C for 2 mins
135 during sampling prior to being heated to 200 °C and backflushed with carrier gas during the desorption cycle. Separation was
136 achieved using a Rtx-624 capillary column (Restek; 105m × 0.32 mm, df = 1.8 µm). The carrier gas was hydrogen and the
137 flow rate was held constant at 1 mL min⁻¹, with a temperature program starting at 30 °C held for 5 minutes, ramped at 25 °C
138 min⁻¹ to 200 °C with a final hold of 25 minutes. The total run time was 42 minutes.

139

140 Figure 1 shows the chromatograms obtained from both instruments (cryo-GC-FID, blue; GC-MS/FID, red) and demonstrates
141 that all compounds, with the exception of 1,2,4-trimethyl benzene (1,2,4-TMB) and 3-carene, and acetone and dimethyl sulfide
142 (DMS), were baseline separated. The chromatogram in Fig. 1 shows a valley between the 1,2,4-TMB and 3-carene peaks and
143 between the acetone and DMS peaks that provides sufficient separation to obtain robust and repeatable peak areas for all four
144 compounds.

145 **Figure 1.** Chromatogram of PRM 0819 showing separation of compounds in the GC-FID/MS (red trace) and cryo-GC-FID
146 (blue trace).
147



148

149 2.3 Stability assessment

150 Stability of all six PRMs were assessed by tracking the ratios of the FID responses of each component relative to an internal
151 reference that was present in every mixture and which is known to be stable (Rhoderick, 2010; Rhoderick and Lin, 2013;
152 Worton et al., 2022). Propane was originally included as an internal tracer to monitor stability but as the PTR-MS in H_3O^+
153 mode cannot detect this compound it was replaced by benzene. Benzene is a good internal tracer with stability of > 2 years
154 that has been well demonstrated relative to propane and hexane for this cylinder type at $5 \mu\text{mol mol}^{-1}$ with an uncertainty of
155 0.5 % (Rhoderick et al., 2019). A similar performance would be expected at $1 \mu\text{mol mol}^{-1}$ and is demonstrated in this work
156 albeit with an uncertainty of 1 % (Fig. S1, Supporting Information).

157

158 All the measurements used for the stability analysis were collected on the same GC-FID/MS instrument with the exception of
159 methanol and acetaldehyde (Cryo-GC-FID). Both instruments remained unchanged throughout the entire time-period of the
160 measurements, which spanned more than 4 years. The observed responses for each compound were corrected for differences
161 in the gravimetric amount fraction and ratioed against the response of the internal reference compound benzene, that was
162 present in every mixture. The uncertainties in the observed ratios included uncertainties for the gravimetric preparation and

163 the repeatability of the analyses. The combined standard uncertainties were multiplied by a coverage factor of 2 ($k=2$) providing
164 a coverage probability of 95 %. The observed ratios were normalised to the average response of all data for that compound to
165 enable comparisons between compounds with different FID responses. For this analysis the data for all six PRMs were
166 considered together to enable an understanding of stability across a longer time period than would be possible for a single
167 PRM. Least squares fit straight-line regressions were modelled to the temporal evolution of the data to determine if there was
168 any statistically significant change in amount fraction of any of the compounds in the PRMs. The slopes from these regression
169 analyses were evaluated with an analysis of variance (ANOVA) test using the 'StatsLinearRegression' function in IGOR pro
170 8.04 (Wavemetrics) (Zar, 1999; Snedecor and Cochran, 1989) to determine whether they were significantly different to zero
171 using an F -test, i.e., no drift in amount fraction during the measurement period ($F < F_c$).
172

173 2.4 Validation approach

174 Five PRMs (0917a, 0917b, 1218, 0119, 0819, and 0821) were validated against PRM 0819 with the exception of PFTBA and
175 toluene that were only present in the most recent PRM (0821). PRM 0819 was used as the reference for all the validations
176 because the parents used for the preparation of this mixture were deliberately different from all other mixtures with the explicit
177 goal of enabling the most robust validation. All compounds were analysed on the GC-FID/MS system, with the exception of
178 methanol and acetaldehyde (Cryo-GC-FID). Toluene was validated by comparison against an existing PRM containing BTEX
179 (benzene, toluene, *m*-xylene, *p*-xylene and *o*-xylene) components that was prepared gravimetrically at NPL in 2018 and had
180 been independently validated against other PRMs that were internationally compared as part of NPL's participation in key
181 comparison CCQM-K10.2018 (Cecelski et al., 2022). These BTEX PRMs are known to be stable for more than 5 years and at
182 the time of the comparisons the BTEX PRM was less than 3 years old. PFTBA was validated against the gravimetric data
183 used to make two independent certified reference materials. The majority of the validation work took place between September
184 and December 2020 with one in 2019 and 2022, respectively, and three in 2021 (Table S3, Supporting Information). As such
185 there is an influence of stability on the validation data as the PRMs differed in age at the time of validation.
186

187 Each comparison was conducted by running the PRMs (0917a, 0917b, 1218, 0119, and 0821) against PRM 0819 in a repeating
188 alternating pattern, (AB)_nA where A represents PRM 0819 and B one of the other PRMs (j) and with the number of repeats
189 ranging between 3 and 5 ($n = 3$ to 5). The ratio in response was determined by dividing B by the average response of the A's
190 immediately before and after each analysis of B. The average ratio was calculated for each compound based on the number of
191 repeats along with the associated standard deviation. The assigned analytical value for compound i in PRM j ($x_{a,i,j}$) was
192 calculated by multiplying the average ratio by the gravimetric amount fraction ($x_{s,i}$) of compound i in PRM 0819. The relative
193 difference (Δx) between the assigned analytical value and the gravimetric value of compound i in PRM j was calculated from:
194

$$\Delta x (\%) = \frac{(x_{u,i,j} - x_{s,i})}{x_{s,i}} \times 100$$

Eq. 2

196

197 The uncertainty in the relative difference combined the standard uncertainty in the repeatability in the analysis with the
198 gravimetric uncertainty. The combined standard uncertainty was multiplied by a coverage factor of 2 ($k=2$) providing a
199 coverage probability of 95 %.

200 **3 Results**

201 **3.1 Composition**

202 The PTR-MS transmission curve reference material contains 20 different VOCs that span a wide range of molecular masses,
203 boiling points, vapour pressures (Table S2, Supporting Information) and functional group classes including alcohol, aldehyde,
204 ketone, alkene, aromatic, halocarbon and siloxane (Table 2). With the PTR-MS technique, most VOCs are entirely detected at
205 their protonated mass, as well as a few compounds that fragment during protonation (e.g. monoterpenes, siloxanes, and
206 isoprene; see Table 1). The compounds were chosen by considering the needs of the PTR-MS user community to cover the
207 full range of mass-to-charge ratios (m/Q) encountered, their low fragmentation following proton transfer and because many
208 are of relevance in atmospheric measurements, which was the initial intended target end user group. Other compounds were
209 included as a consequence of the preparation method, that is the case for *n*-hexane, which is present as the solvent for D3-
210 siloxane, and propane, which was present in one of the parent mixtures and was originally included as an internal tracer to
211 monitor stability. The composition evolved over time, as shown in Table 2, with DMS, 1,2,4-TCB, D3-siloxane, toluene and
212 PFTBA being added at different times, and propane being removed in the final iteration. For D4-siloxane there was a
213 preparation error, and it was not added to either PRM 1218 or 0119.

214

215 An amount fraction of nominally 1 $\mu\text{mol mol}^{-1}$ in a balance gas of nitrogen was selected as a compromise between preparation
216 complexity and mixture stability. This amount fraction enabled many components to be prepared from parent mixtures of
217 higher amount fraction ($\geq 10 \mu\text{mol mol}^{-1}$), which substantially simplifies the preparation process. This amount fraction also
218 provided a reasonable starting point for stability of the wide range of function groups present in the mixture some of which are
219 known to have more limited stability at lower abundances fractions (nmol mol^{-1}) (Allen et al., 2018).

220 **3.2 Traceability to the International System of Units (SI)**

221 Traceability of the primary realisations to the international community through CCQM key comparisons or regional
222 EURAMET comparisons provides confidence in the accuracy of the amount fractions for all components. SI traceability is
223 important for underpinning long-term measurements as it provides a stable anchor point with which to reference all
224 measurements to. Table 1 shows which comparisons underpin the traceability for each of the different components. All the

225 components are underpinned by at least one CCQM or EURAMET comparison with the exception of 1,2,4-TFB, 1,2,4-TCB
226 and PFTBA, for which there are currently no existing relevant comparisons.

227 3.3 Hierarchies

228 Table 2 shows all the parent mixtures and their preparation dates used to prepare all six PRMs (0917a, 0917b, 1218, 0119,
229 0819, and 0821) and in total 50 different parent mixtures were used. In general, parent mixtures were similar for PRMs 0917a,
230 0917b, 1218 and 0119 but were different to PRMs 0819 and 0821 providing independence and thus confidence in the validation
231 work and in the preparations. There were a few exceptions. For *m*-xylene the parent used for PRMs 1218 and 0119 was the
232 same as PRM 0821 but different to 0917a, 0917b and 0821. For 1,2,4-TMB only two parent mixtures were used one for 0917a,
233 0917b and 0821 and another for 1218, 0119 and 0819. For 3-carene only two parents were used one for 0917a and 0917b and
234 another for 1218, 0119, 0819 and 0821. For D3-siloxane three parents were used, one for 1218 and 0819, one for 0119 and
235 another for 0821.

236 3.4 PRM Validation

237 Figure 2 shows the relative differences (Δx) determined from Eq. 1 for all compounds using all the validation data obtained
238 from the 13 comparisons outlined in Table S3 (Supporting Information). In the majority of cases PRM 0819 was used as the
239 reference to which all others are compared. It was chosen as such because at that time it was the newest PRM to be produced
240 and was used to benchmark all the others that had already been made. Thus, PRM 0821 was also referenced to PRM 0819 to
241 provide a link between all six PRMs. All the data shown in Fig. 2 is the FID data from the GC-MS/FID instrument with the
242 exception of acetonitrile (MS data from the GC-MS/FID instrument), methanol and acetaldehyde (FID data from the cryo-GC-
243 FID instrument). The MS data is used for acetonitrile because the FID data shows a larger variability, which is likely attributed
244 to the co-elution of an impurity in the FID that was present at different amount fractions in the different PRMs but we do not
245 have an conclusive evidence to support this and additional work is needed to confirm. This variability is not observed in the
246 MS data providing better precision (Fig. S2, Supporting Information).

247
248 In general, the data from Fig. 2 could be split into three groups. The first group consisted of propane, isoprene, benzene,
249 toluene, 3-carene, methanol, acetonitrile, acetaldehyde, *m*-xylene, 1,2,4-TMB, MEK where the spread in the validation data is
250 within 3 % and these represent components where NPL had substantial prior experience. The second group is acetone, DMS,
251 MVK and PFTBA where the spread in the validation data is within 5 % and these are relatively new components where
252 capabilities were developed more recently. Recognising the challenges in preparing PRMs containing siloxanes as a result of
253 their lower vapour pressures and observing the recent improvements in preparation since 2019, the D4-siloxane and D5-
254 siloxanes can also be categorised as group 2 after excluding the earliest parent preparations used for 0917a and 0917b in 2017,
255 which are inconsistent with more recent work as part of the EURAMET 1305 Siloxanes comparison (Van Der Veen et al.,
256 2022). The final group is comprised of D3-siloxane and 1,2,4-TCB where the spread in validation data is within 10 % and

257 these compounds represent those which the most challenging to prepare as a result of either unique phase transition properties
258 or low vapour pressures, respectively. There is an observable bias of about 8 % between two groups of mixtures; one group is
259 1218 and 0819 and the other is 0119 and 0821. This reflects differences between the parent mixtures (2586, 2693 and 3134)
260 that resulted from the challenges in preparation. Ethanol also sits with this group in part due to the small size of the peak
261 observed in the GC-MS/FID instrument and because of what looks like an outlier (0119), suggesting some potential losses
262 during preparation that were unique to this one PRM.

263

264 All the FID and supporting MS data for all compounds are shown in Fig. S2 (Supporting Information). No MS data was
265 available for toluene, 1,2,4-TCB or PFTBA because the relevant single m/Q ions had not been included in the MS single ion
266 monitoring method at the time of analysis and methanol where the MS signal was too small to provide a reliable response.
267 Figure S2 shows very good agreement between the FID and MS validation with all components agreeing within the
268 uncertainties providing confidence in the validation results.

269


270 In addition to the observed bias in parent mixtures for D3-siloxane three other parent mixtures were also discovered to be
271 biased after re-analysis. The observed differences have been corrected for in Fig. 2 and Fig. S2 (Supporting Information). For
272 methanol, one parent (A410) was confirmed to be 5.0 % high relative to the other parents (A463, A540 and A602) For MVK,
273 one parent (3070) was confirmed to be 6.3 % low relative to the other parents (2064 and 2088). For 1,2,4-TMB, one parent
274 D711530 was confirmed to be 6.0 % low relative to D442684 and other in-house standards of 1,2,4-TMB not used in this work
275 but used to prepare 30 component ozone precursor mixtures at NPL (Grenfell et al., 2010).

276

277

278 **Figure 2.** Relative difference (Δx) using the FID data (except acetonitrile, which uses the MS data) for all components in five
 279 of the PRMs (0917a, 0917b, 1218, 0119 and 0821) relative to PRM 0819 (solid symbols). The solid black line represents the
 280 average of these validations with the error bar representing the associated expanded uncertainty (2σ). For D4- and D5-siloxane
 281 the averages do not include the validations from 0917a or 0917b. Methanol and acetaldehyde data are from the cryo-GC-FID
 282 instrument while all others are from the GC-FID/MS instrument. Open symbols represent the original data before correcting
 283 for biases observed in three of the parent mixtures (A410, 5 % low for methanol; 3070, 6.3 % low for MVK and D711530 6
 284 % low for 1,2,4-TMB). PFTBA and toluene were only included in the most recently prepared PRM (0821) and are not present
 285 in 0819. Their validation is described in the text. Supporting validation data from all the MS and FID measurements is shown
 286 in Fig. S2 (Supporting Information).

287



288

289

290 3.5 CRM Validation

291 To enable a more cost effective and timely delivery to end users a certified reference material (CRM) was also developed. In
 292 contrast to the PRMs the CRMs are not prepared by gravimetry but by the direct addition of multicomponent mixtures derived
 293 from the original pure liquids. Further details of the preparation method are given in the Supporting Information
 294 (Supplementary text: Preparation and validation of certified reference materials). The amount fractions for the components in

295 the CRMs were assigned through analytical comparisons between each CRM and one or more of the PRMs. In this way,
296 preparation is quicker and more cost effective while maintaining the integrity of the values and their traceability. An additional
297 advantage of the CRMs is that because the solid D3-siloxane is dissolved in the other components no *n*-hexane is used which
298 avoids any potential interferences from the presence of reagent ions other than H_3O^+ like O_2^+ and NO^+ . Initially with the
299 developed CRM preparation method it was possible to produce mixtures that had blend tolerances of 20 – 30 % (Fig. S3,
300 Supporting Information), which are suitable for end users but work is continuing to improve this with the aim of achieving
301 better than 10 % blend tolerances in the near future. The blend tolerances are just an indication of the repeatability of the
302 preparation process and do not reflect the uncertainties in the assigned value, which are between 3 – 10 % (compound
303 dependent). These uncertainties were dominated by the observed differences between the PRMs.

304 3.6 Stability

305 Figure 3 shows stability data for four selected compounds; methanol, isoprene, D3-siloxane and PFTBA. These were selected
306 as representative examples of the different observed stability behaviours although the stability data plots corresponding to all
307 compounds are shown in Fig. S1 (Supporting Information). The trendlines from the least squares fit straight-line regressions
308 shown in Fig. 3 and Fig. S1 (Supporting Information) were used to determine the annual drift rates shown in Table 3 and Fig.
309 4. From the ANOVA test there are statistically significant trends ($F > F_c$) for 10 of the compounds (methanol, acetonitrile,
310 acetaldehyde, ethanol, acetone, DMS, isoprene, MVK, benzene, D3-siloxane) but these trends are small ($< 2\% \text{ yr}^{-1}$) except
311 for methanol, acetonitrile and PFTBA.

312


313 Methanol and acetaldehyde were the only two components that were measured on the cryo-GC-FID and hence the datasets are
314 more limited. A result is that there is no overlap between the three PRMs so any systematic differences between them may
315 result in an artificial bias, which may exacerbate any stability trend. More work is needed to confirm this. The stability data
316 for D3-siloxane reflects the observed validation bias and shows two clear trends; one for 1218 and 0819 and the other for 0119
317 and 0821. For the regression analysis and drift calculations these have been treated independently (Table 3).

318

319 All compounds, with the exceptions of methanol, acetonitrile and PFTBA, show trends similar to isoprene with good stability
320 and annual drift rates of $< 3\% \text{ yr}^{-1}$ (Table 3). For acetonitrile the large spread in validation data (FID data; Fig. S2, Supporting
321 Information) leads to a noisy stability dataset that may play a role in the larger observed drift rate or this component maybe
322 less stable. As PFTBA was only included in the last PRM (0821) the stability data only represents about half a year and
323 extrapolating the current trend to 1 and 2 years results in a drift rate that is not accurate as interpolation of the data would
324 suggest no statistical change in amount fraction and minimal drift. More data is needed to confirm the longer-term stability
325 behaviour of PFTBA.

326

327 **Figure 3.** Stability of normalised response with time for four selected compounds relative to benzene, methanol (top left),
 328 isoprene (top right), D3-siloxane (bottom left) and PFTBA (bottom right) for all six PRMs (solid symbols). The open symbols
 329 (methanol; top left) show the original data before being corrected for an observed 5.0 % bias in the parent mixture (A410).
 330 The best fit curves from least squares straight line regression analyses are shown (solid black line) along with the 95 %
 331 confidence interval of the fits (shaded area). The slope, intercept and F-statistic data from the regression analyses are shown
 332 in Table 3. Stability plots for all compounds are shown in Fig. S1 (Supporting Information).
 333



334
 335

336 Given the age differences between the different PRMs at the time of validation (233 – 709 days; Table S2, Supporting
 337 Information) it is not possible to deconvolute the contributions of stability and preparation to the observed validation
 338 differences. However, Fig. 4 shows that for the majority of compounds there is good agreement between the observed
 339 average validation data and the calculated drift for over 1 – 2 years, with the exception of methanol, acetonitrile and PFTBA,
 340 which differ for the reasons discussed previously. These observations are consistent with the age differences of the different
 341 PRMs at the time of validation indicating that stability was likely the major driver between the observed validation
 342 differences.

343 **Table 3.** Summary of the results of the least squares straight-line regression analysis for all stability data shown in Fig. 3 and
 344 Fig. S2 (Supporting Information). Results are shown for the slope ($\pm 2\sigma$), intercept ($\pm 2\sigma$), ANOVA test statistics (F and Fc)
 345 used to evaluate the presence of a statistically significant trend ($F > F_c$), the calculated annual drift ($\pm 2\sigma$) determined from the
 346 linear fit and the average of the validation data ($\pm 2\sigma$), also shown in Fig. 2.

347


compound	slope ($\times 10^{-5}$)	intercept	F	Fc	Drift (%/yr)	Avg valid. (%)
methanol ^a	-9.539 ± 2.700	1.067 ± 0.021	57.005	5.318	-3.48 ± 0.42	-0.37 ± 2.77
acetonitrile	-12.328 ± 6.128	1.036 ± 0.026	16.530	4.085	-4.50 ± 0.94	-0.72 ± 2.08
acetaldehyde ^a	-5.345 ± 2.800	1.037 ± 0.022	19.699	5.318	-1.95 ± 0.44	-0.40 ± 2.53
propane	0.653 ± 5.393	0.997 ± 0.027	0.062	4.225	0.24 ± 1.97	0.16 ± 1.66
ethanol	-7.841 ± 7.55	1.023 ± 0.032	4.405	4.085	-2.86 ± 0.36	-0.61 ± 5.64
acetone	3.462 ± 3.206	0.990 ± 0.013	4.765	4.085	1.26 ± 0.86	1.18 ± 3.08
DMS	2.441 ± 2.351	0.995 ± 0.007	4.473	4.149	0.89 ± 2.24	-0.76 ± 3.22
isoprene	-1.338 ± 0.975	1.004 ± 0.004	7.690	4.085	-0.49 ± 1.17	-0.04 ± 0.60
MVK	-3.523 ± 2.564	1.010 ± 0.011	7.708	4.085	-1.29 ± 0.94	-0.61 ± 4.50
MEK	0.575 ± 1.967	0.998 ± 0.008	0.349	4.085	0.21 ± 0.36	-1.23 ± 2.25
Benzene ^b	1.329 ± 0.983	0.996 ± 0.004	7.456	4.085	0.49 ± 0.18	-0.48 ± 1.07
Toluene ^c	-3.546 ± 4.536	1.002 ± 0.004	2.902	4.747	-1.30 ± 1.66	0.19 ± 0.29
m-xylene	0.129 ± 2.034	1.000 ± 0.009	0.016	4.085	0.05 ± 0.74	-0.87 ± 1.88
1,2,4-TMB	-0.870 ± 5.155	1.003 ± 0.022	0.116	4.085	-0.32 ± 1.69	-0.57 ± 2.42
1,2,4-TFB	-1.373 ± 1.448	1.004 ± 0.006	3.672	4.085	-0.50 ± 2.05	0.27 ± 0.74
+3-carene	-0.734 ± 4.631	1.002 ± 0.019	0.103	4.085	-0.27 ± 2.84	-0.25 ± 1.33
1,2,4-TCB	4.512 ± 6.455	0.991 ± 0.018	2.027	4.149	1.65 ± 1.16	-1.73 ± 6.56
D3-siloxane ^d	-2.641 ± 1.740	1.056 ± 0.007	11.444	4.965	-0.96 ± 0.29	-4.02 ± 9.67
D4-siloxane ^e	3.195 ± 3.220	0.970 ± 0.007	4.287	4.351	1.17 ± 0.56	-4.02 ± 9.67
D5-siloxane ^e	4.799 ± 4.300	0.988 ± 0.012	0.765	4.225	1.75 ± 0.74	-2.03 ± 4.06
D5-siloxane ^e	2.066 ± 0.390	0.985 ± 0.026	1.833	4.085	0.75 ± 1.68	-0.49 ± 3.27
PFTBA ^c	-12.045 ± 13.440	1.007 ± 0.010	3.813	4.747	-4.40 ± 1.66	3.31 ± 0.70

348 ^aThe GC-FID data for methanol and acetaldehyde was too small to be quantified so this data is from the cryo-GC-FID data and is limited. ^bBenzene stability
 349 was determined relative to isoprene. ^cToluene and PFTBA were only included in the most recent PRM so the assessment of stability is limited in its duration
 350 to only 200 days. ^dThere was a clear bias between several of the PRMs caused by differences in the parent mixtures used so the trends were fitted to the two
 351 obvious groupings. ^eData from 0917a and 0917b were excluded from the regression analysis.

352

353 **Figure 4.** Comparison of 1-year (filled grey squares) and 2-year (open grey squares) drift rates, calculated from the data in
354 Table 3, with the average validation data (black bars) taken from Fig. 2. For D3-siloxane there are two datapoints for the drift
355 correspond to the two regressions shown in Table 3. The error bars represent the associated expanded uncertainties,
356 representing the 95 % confidence limit.

357



358

359 4 Conclusions

360 In this work the development of new primary reference materials (PRMs) and certified reference materials (CRMs) for
361 constraining the mass dependent transmission curve of PTR-MS instruments have been described along with an evaluation of
362 the validation and stability of the PRMs and the repeatability in preparation (blend tolerances) for the CRMs. Six of these
363 PRMs have been prepared to date from a suite of 50 parent mixtures and these have been used to value assign more than 10
364 CRMs that have been disseminated to end users. In general, there is evidence of very good agreement for the majority of
365 components that supports the robustness of the preparation and 2 years of stability. Challenges were observed in preparation
366 for the least volatile compounds especially for D3-siloxane due to it being a solid at room temperature and pressure. More
367 work is needed to better describe the long-term stability of methanol, acetonitrile and PFTBA. This work highlighted several
368 challenges in analysis that could be resolved by the development of a new analytical method utilised a single instrument
369 equipped with both a preconcentration trap and dual detector setup (MS and FID). This work demonstrates what is currently

370 possible with respect to composition, amount fraction, uncertainty and stability and provides an important reference to which
371 other gas standards that are in use with the PTR-MS can be compared and benchmarked to verify their accuracy to further
372 improve the comparability of PTR-MS measurement data.

373

374 In the short term (next 5 years) the implementation of an SI traceable transmission curve reference material, such as the one
375 described in this work, using a method similar to that described in Holzinger et al., 2019 is the most pragmatic approach to
376 directly address improving the accuracy of quantitation and comparability between different PTR-MS instruments and users.
377 This reflects the challenges and complications of rapidly developing a universally accepted calibration system based on pure
378 liquids that is SI traceable. The use of a SI traceable reference material to properly constrain the transmission curve provides
379 a readily applicable framework to ensure confidence in temporal and spatial data to support the use of PTR-MS in a broad
380 range of application areas. The use of the transmission curve reference material approach should be seen as a pre-requisite and
381 a complement to additional future efforts to provide alternative calibration efforts for specific target compounds where
382 uncertainties of better than 30 % are needed. Alternative approaches would certainly be necessary for those compounds that
383 are unsuitable for inclusion in high pressure gas standards possibly as a result of very low vapour pressures or other
384 complicating factors such as chemical compatibility with other compounds.

385

386 Future work to improve the uncertainty of individual components that have the greatest influence on the transmission curve fit
387 would have the biggest influence on the accuracy and repeatability of the transmission curve retrieval thereby maximising the
388 impact of future improvements for the PTR-MS user community. For PTR-MS instruments that utilise time of flight mass
389 spectrometers the focus would be on improving the uncertainty of the largest molecular weight components specifically the
390 D3-, D4-, D5-siloxanes and 1,2,4-TCB, which represent the greatest challenges in preparation due to their low vapour
391 pressures.

392 **Data availability**

393 All data used to produce the figures in this paper are available on request.

394 **Competing interests**

395 The authors declare that they have no conflict of interest.

396 **Author contributions**

397 DRW and RH conceptualized the work. DRW processed the data, produced the figures and tables and wrote the paper. RH
398 provided inputs to define the composition of the PTR-MS reference material and contributed to the writing of the paper. SM
399 developed novel methods for the preparation of primary and certified reference materials, planned and prepared all reference
400 materials, conducted all the validation and stability data collection and contributed to the writing of the paper. KOD contributed
401 to the preparation of reference materials, reprocessed some of the stability data and worked with SM to prepare the certified
402 reference materials. All authors reviewed the paper.

403 **Financial support**

404 NPL acknowledges funding from the Department of Business Energy and Industrial Strategy (BEIS) National Measurement
405 System. RH acknowledges funding from the EMPIR programme (19ENV06 MetClimVOC), co-financed by the Participating
406 States and from the European Union's Horizon 2020 research and innovation programme, and from the European Union's
407 Horizon 2020 research and innovation programme (ACTRIS-2) under grant agreement no. 654109.

408 **References**

- 409 Allen, N. D. C., Worton, D. R., Brewer, P. J., Pascale, C., and Niederhauser, B.: The importance of cylinder passivation for preparation and
410 long-term stability of multicomponent monoterpene primary reference materials, *Atmos. Meas. Tech.*, 11, 6429-6438, 10.5194/amt-11-6429-
411 2018, 2018.
- 412 Ammann, C., Spirig, C., Neftel, A., Steinbacher, M., Komenda, M., and Schaub, A.: Application of PTR-MS for measurements of biogenic
413 VOC in a deciduous forest, *International Journal of Mass Spectrometry*, 239, 87-101, 10.1016/j.ijms.2004.08.012, 2004.
- 414 Beauchamp, J., Herbig, J., Dunkl, J., Singer, W., and Hansel, A.: On the performance of proton-transfer-reaction mass spectrometry for
415 breath-relevant gas matrices, *Measurement Science and Technology*, 24, 125003, 2013.
- 416 Biasioli, F., Gasperi, F., Yeretzyan, C., and Märk, T. D.: PTR-MS monitoring of VOCs and BVOCs in food science and technology, *TrAC*
417 *Trends in Analytical Chemistry*, 30, 968-977, <https://doi.org/10.1016/j.trac.2011.03.009>, 2011.
- 418 Blake, R. S., Monks, P. S., and Ellis, A. M.: Proton-Transfer Reaction Mass Spectrometry, *Chemical Reviews*, 109, 861-896,
419 10.1021/cr800364q, 2009.
- 420 Brown, A. S., Milton, M. J. T., Brookes, C., Vargha, G. M., Downey, M. L., Uehara, S., Augusto, C. R., Fioravante, A. d. L., Sobrinho, D.
421 G., Dias, F., Woo, J. C., Kim, B. M., Kim, J. S., Mace, T., Fükö, J. T., Qiao, H., Guenther, F., Rhoderick, J., Gameson, L., Botha, A.,
422 Tshilongo, J., Ntsasa, N. G., Val'ková, M., Durisova, Z., Kustikov, Y., Konopelko, L., Fatina, O., and Wessel, R.: Final report on CCQM-
423 K93: Preparative comparison of ethanol in nitrogen, *Metrologia*, 50, 08025-08025, 10.1088/0026-1394/50/1a/08025, 2013.
- 424 Brown, R. J. C. and Milton, M. J. T.: Developments in accurate and traceable chemical measurements, *Chemical Society Reviews*, 36, 904-
425 913, 10.1039/b507452p, 2007.
- 426 Cappellin, L., Karl, T., Probst, M., Ismailova, O., Winkler, P. M., Soukoulis, C., Aprea, E., Märk, T. D., Gasperi, F., and Biasioli, F.: On
427 Quantitative Determination of Volatile Organic Compound Concentrations Using Proton Transfer Reaction Time-of-Flight Mass
428 Spectrometry, *Environmental Science & Technology*, 46, 2283-2290, 10.1021/es203985t, 2012.
- 429 Cecelski, C. E., Rhoderick, G. C., Possolo, A. M., Carney, J., Vokoun, M., Privoznikova, J., Lee, S., Kang, J. H., Kim, Y. D., Kim, D. H.,
430 Macé, T., Sutour, C., Pascale, C., Ntsasa, N., Tshilongo, J., Jozela, M., Leshabane, N., Lekoto, G., Worton, D. R., Brewer, P. J., Farrow-
431 Dunn, F., Moreno, S., Wirtz, K., Stummer, V., Konopelko, L. A., Kolobova, A. V., Kustikov, Y. A., Klimov, A. Y., Efremova, O. V., Wijk,
432 J. I. T. v., and Veen, A. M. H. v. d.: International comparison CCQM-K10.2018: BTEX in nitrogen at 5 nmol mol⁻¹,
433 *Metrologia*, 59, 08003, 10.1088/0026-1394/59/1a/08003, 2022.
- 434 de Gouw, J. and Warneke, C.: Measurements of volatile organic compounds in the earth's atmosphere using proton-transfer-reaction mass
435 spectrometry, *Mass Spectrom. Rev.*, 26, 223-257, 10.1002/mas.20119, 2007.

436 de Gouw, J., Warneke, C., Karl, T., Eerdeken, G., van der Veen, C., and Fall, R.: Sensitivity and specificity of atmospheric trace gas
437 detection by proton-transfer-reaction mass spectrometry, *International Journal of Mass Spectrometry*, 223, 365-382, 10.1016/s1387-
438 3806(02)00926-0, 2003.

439 Grenfell, R. J. P., Brookes, C., Vargha, G., Quincey, P., Milton, M., Woods, P. T., and Harris, P.: Euramet 886 Comparison of multi-
440 component ambient VOC measurements, National Physical Laboratory, 2008.

441 Grenfell, R. J. P., Milton, M. J. T., Harling, A. M., Vargha, G. M., Brookes, C., Quincey, P. G., and Woods, P. T.: Standard mixtures of
442 ambient volatile organic compounds in synthetic and whole air with stable reference values, *Journal of Geophysical Research-Atmospheres*,
443 115, 10.1029/2009jd012933, 2010.

444 Hansel, A., Jordan, A., Warneke, C., Holzinger, R., Wisthaler, A., and Lindinger, W.: Proton-transfer-reaction mass spectrometry (PTR-
445 MS): on-line monitoring of volatile organic compounds at volume mixing ratios of a few pptv, *Plasma Sources Science & Technology*, 8,
446 332-336, 10.1088/0963-0252/8/2/314, 1999.

447 Holzinger, R., Lee, A., McKay, M., and Goldstein, A. H.: Seasonal variability of monoterpene emission factors for a ponderosa pine
448 plantation in California, *Atmos. Chem. Phys.*, 6, 1267-1274, 10.5194/acp-6-1267-2006, 2006.

449 Holzinger, R., Acton, W. J. F., Bloss, W. J., Breitenlechner, M., Crilley, L. R., Dusanter, S., Gonin, M., Gros, V., Keutsch, F. N., Kiendler-
450 Scharr, A., Kramer, L. J., Krechmer, J. E., Languille, B., Locoge, N., Lopez-Hilfiker, F., Materić, D., Moreno, S., Nemitz, E., Quéléver, L.
451 L. J., Sarda Esteve, R., Sauvage, S., Schallhart, S., Sommariva, R., Tillmann, R., Wedel, S., Worton, D. R., Xu, K., and Zaytsev, A.: Validity
452 and limitations of simple reaction kinetics to calculate concentrations of organic compounds from ion counts in PTR-MS, *Atmos. Meas.*
453 *Tech.*, 12, 6193-6208, 10.5194/amt-12-6193-2019, 2019.

454 ISO: 6142-1:2015 Gas analysis - Preparation of calibration gas mixtures - Gravimetric method for class I mixtures. International Organization
455 for Standardisation (ISO), Geneva, Switzerland, 2015.

456 ISO: 19229:2019 Gas analysis - Purity analysis and the treatment of purity data. International Organization for Standardisation (ISO),
457 Geneva, Switzerland, 2019.

458 Lee, S., Kang, J. H., Kim, Y. D., Kim, D. H., Jung, J., Ma, H., Wu, H., Bartlett, S., Worton, D., Murugan, A., Brewer, P. J., Konopelko, L.
459 A., Kolobova, A. V., Malginov, A. V., Dobryakov, Y. G., Pankratov, V. V., and Efremova, O. V.: International comparison CCQM-K165:
460 dimethyl sulfide in nitrogen at 5 nmol mol⁻¹, *Metrologia*, 59, 08011, 10.1088/0026-1394/59/1A/08011, 2022.

461 Lee, S., Heo, G. S., Kim, Y., Oh, S., Han, Q., Wu, H., Konopelko, L. A., Kustikov, Y. A., Kolobova, A. V., Efremova, O. V., Pankratov, V.
462 V., Pavlov, M. V., Culleton, L. P., Brown, A. S., Brookes, C., Li, J., Ziel, P. R., and van der Veen, A. M. H.: International key comparison
463 CCQM-K94: 10 μmol/mol dimethyl sulfide in nitrogen, *Metrologia*, 53, 08002, 10.1088/0026-1394/53/1a/08002, 2016.

464 Liaskos, C., Rhoderick, G., Hodges, J., Possolo, A., Goodman, C., Kim, Y. D., Kim, D. H., Lee, S., Allen, N., Corbel, M., Worton, D.,
465 Brown, R., and Brewer, P.: CCQM-K121 - Monoterpenes in nitrogen at 2.5 nmol mol⁻¹ final report, *Metrologia*, 55, 08019,
466 10.1088/0026-1394/55/1a/08019, 2018.

467 Lindinger, W., Hansel, A., and Jordan, A.: On-line monitoring of volatile organic compounds at pptv levels by means of proton-transfer-
468 reaction mass spectrometry (PTR-MS) medical applications, food control and environmental research, *International Journal of Mass*
469 *Spectrometry and Ion Processes*, 173, 191-241, [http://dx.doi.org/10.1016/S0168-1176\(97\)00281-4](http://dx.doi.org/10.1016/S0168-1176(97)00281-4), 1998.

470 Müller, M., Mikoviny, T., and Wisthaler, A.: Detector aging induced mass discrimination and non-linearity effects in PTR-ToF-MS,
471 *International Journal of Mass Spectrometry*, 365-366, 93-97, <https://doi.org/10.1016/j.ijms.2013.12.008>, 2014.

472 Park, J. H., Goldstein, A. H., Timkovsky, J., Fares, S., Weber, R., Karlik, J., and Holzinger, R.: Eddy covariance emission and deposition
473 flux measurements using proton transfer reaction – time of flight – mass spectrometry (PTR-TOF-MS): comparison with PTR-
474 MS measured vertical gradients and fluxes, *Atmos. Chem. Phys.*, 13, 1439-1456, 10.5194/acp-13-1439-2013, 2013.

475 Rhoderick, G. C.: Stability assessment of gas mixtures containing terpenes at nominal 5 nmol/mol contained in treated aluminum gas
476 cylinders, *Analytical and Bioanalytical Chemistry*, 398, 1417-1425, 10.1007/s00216-010-4058-0, 2010.

477 Rhoderick, G. C. and Lin, J.: Stability Assessment of Gas Mixtures Containing Monoterpenes in Varying Cylinder Materials and Treatments,
478 *Analytical Chemistry*, 85, 4675-4685, 10.1021/ac400324v, 2013.

479 Rhoderick, G. C., Cecelski, C. E., Miller, W. R., Worton, D. R., Moreno, S., Brewer, P. J., Viallon, J., Idrees, F., Moussay, P., Kim, Y. D.,
480 Kim, D., Lee, S., Baldan, A., and Li, J.: Stability of gaseous volatile organic compounds contained in gas cylinders with different internal
481 wall treatments, *Elementa-Sci. Anthropol.*, 7, 22, 10.1525/elementa.366, 2019.

482 Snedecor, G. W. and Cochran, W. G.: *Statistical methods*, Iowa State University Press 1989.

483 Steinbacher, M., Dommen, J., Ammann, C., Spirig, C., Neftel, A., and Prevot, A. S. H.: Performance characteristics of a proton-transfer-
484 reaction mass spectrometer (PTR-MS) derived from laboratory and field measurements, *International Journal of Mass Spectrometry*, 239,
485 117-128, <https://doi.org/10.1016/j.ijms.2004.07.015>, 2004.

486 Taipale, R., Ruuskanen, T. M., Rinne, J., Kajos, M. K., Hakola, H., Pohja, T., and Kulmala, M.: Technical Note: Quantitative long-term
487 measurements of VOC concentrations by PTR-MS – measurement, calibration, and volume mixing ratio calculation methods, *Atmos.*
488 *Chem. Phys.*, 8, 6681-6698, 10.5194/acp-8-6681-2008, 2008.

489 Van der Veen, A. M. H., de Krom, I., Nieuwenkamp, G., Culleton, L. P., Worton, D. R., and Allden, J.: Bilateral comparison VSL-NPL on
490 primary gas standards for the siloxane content in biomethane, TBD, 2022.

491 Veen, A. M. H. v. d., Hout, J. W. v. d., Ziel, P. R., Oudwater, R. J., Fioravante, A. L., Augusto, C. R., Brum, M. C., Uehara, S., Akima, D.,
492 Bae, H. K., Kang, N., Woo, J.-C., Liaskos, C. E., Rhoderick, G. C., Jozela, M., Tshilongo, J., Ntsasa, N. G., Botha, A., Brewer, P. J., Brown,
493 A. S., Bartlett, S., Downey, M. L., Konopelko, L. A., Kolobova, A. V., Pankov, A. A., Orshanskaya, A. A., and Efremova, O. V.: International
494 comparison CCQM-K111—propane in nitrogen, *Metrologia*, 54, 08009, 10.1088/0026-1394/54/1a/08009, 2017.
495 Worton, D. R., Moreno, S., Brewer, P. J., Li, J., Baldan, A., and van der Veen, A. M. H.: Bilateral comparison of primary reference materials
496 (PRMs) containing methanol, ethanol and acetone in nitrogen, *Accreditation and Quality Assurance*, 10.1007/s00769-022-01513-y, 2022.
497 Yuan, B., Koss, A. R., Warneke, C., Coggon, M., Sekimoto, K., and de Gouw, J. A.: Proton-Transfer-Reaction Mass Spectrometry:
498 Applications in Atmospheric Sciences, *Chemical Reviews*, 117, 13187-13229, 10.1021/acs.chemrev.7b00325, 2017.
499 Zar, J. H.: *Biostatistical Analysis*, 4th, Prentice Hall, Englewood Cliffs, New Jersey 1999.
500 Zhao, J. and Zhang, R.: Proton transfer reaction rate constants between hydronium ion (H₃O⁺) and volatile organic compounds, *Atmospheric*
501 *Environment*, 38, 2177-2185, <http://dx.doi.org/10.1016/j.atmosenv.2004.01.019>, 2004.
502