

Version 1

EPIC channel	Calibration coefficients published on 16 July 2016	M/E ratio / regression difference (%)
443 nm	8.80×10^{-6}	2.79
551 nm	6.90×10^{-6}	1.98
680 nm	1.00×10^{-5}	1.01
780 nm	1.50×10^{-5}	0.41

Version 2

EPIC channel	Calibration coefficients published on 6 July 2017	M/E ratio / regression difference (%)
443 nm	8.34×10^{-6}	0.1
551 nm	6.66×10^{-6}	0.5
680 nm	9.30×10^{-6}	0.5
780 nm	1.435×10^{-5}	1.4